



# Polypropylene Storage Plates

Microplates

Deep-Well Plates

KingFisher™ Plates

Sample Library Tubes Strip

Sealing Mat




Designed for Safe Storage and High-throughput Screening



# Polypropylene Plates

Microplates | Deep-Well Plates | KingFisher™ Plates | Sample Library Tubes strip | Sealing Mat

- Plates are designed according to ANSI/SLAS format for high throughput screening, sample collection, compound storage, compound libraries, and combinatorial chemistry application.
- Compatible with most multichannel pipettes and automated Liquid-handling equipment.
- Made of virgin polypropylene, conforming USP class VI and compatible with solvents such as DMSO, phenol, and chloroform.
- Manufactured in ISO class 8 cleanroom facilities. Certified free of detectable DNase, RNase, and endotoxin.
- Alphanumeric well coding moulded in A-H and 1-12 axis.
- Raised rims on each well facilitate secure sealing. Compatible with adhesive seals or Silicone Sealing Mats.
- Temperature range for use: -80 °C to 120 °C.



## MICRO PLATES

### 96 Well Standard Volume Micro Plates

Cat #	Well Volume	Height (mm)	Well Shape	Well Bottom	Colour	Packing unit/Case
510060	390µl	14.5	Round	Flat	Natural	5/50
510060BK	390µl	14.5	Round	Flat	Black	5/50
510063	390µl	14.5	Round	Flat	White	5/50
510061	350µl	14.5	Round	U	Natural	5/50
510061BK	350µl	14.5	Round	U	Black	5/50
510064	350µl	14.5	Round	U	White	5/50
510062	340µl	14.5	Round	V	Natural	5/50
510062BK	340µl	14.5	Round	V	Black	5/50
510065	340µl	14.5	Round	V	White	5/50

### 96 Well Extended Volume Micro Plates

520071	500µl	14.5	Round	V	Natural	5/50
520072	500µl	14.5	Round	U	Natural	5/50

### 384 Well Micro Plates

★ 520073	130µl	14.5	Square	U	Natural	5/50
★ 520074	130µl	14.5	Square	Pyramid	Natural	5/50

## DEEP WELL PLATES

### 96 Well Low-profile Deep Well Plate for use with Magnetic Separator

Cat #	Well Volume	Height (mm)	Well Shape	Well Bottom	Colour	Packing unit/ Case
510066	0.8 ml	31	Round	V	Natural	5/50

### 96 Well Low-profile Deep Well Plate

510067	1.1 ml	27	Round	U	Natural	5/50
--------	--------	----	-------	---	---------	------

### 96 Well Standard-profile Deep Well Plates

510069	1.2 ml	41.5	Round	U	Natural	5/50
510068	2.0 ml	45	Round	U	Natural	5/50
510077	2.2 ml	44	Square	Pyramid	Natural	5/50
★ 510078	2.4 ml	44	Square	U	Natural	5/50

### 384 Well Low-profile Deep Well Plate

510220	240µl	22	Square	Pyramid	Natural	5/50
--------	-------	----	--------	---------	---------	------

★ Delivery May'24

# Polypropylene Plates

Microplates | Deep-Well Plates | KingFisher™ Plates | Sample Library Tubes strip | Sealing Mat

## HIGH VOLUME DEEP WELL PLATES

### 48 Well, 5 mL Deep Well Plates

Cat #	Well Volume	Height (mm)	Well Shape	Well Bottom	Colour	Packing unit / Case
510080	5 mL	44	Rectangular	Pyramid	Natural	5/50
510080S	5 mL Sterile	44	Rectangular	Pyramid	Natural	5/50

### 24 Well, 10 mL Deep Well Plate

510081	10 mL	44	Square	Pyramid	Natural	5/50
510081S	10 mL Sterile	44	Square	Pyramid	Natural	5/50
510082	10 mL	44	Square	U	Natural	5/50
510082S	10 mL Sterile	44	Square	U	Natural	5/50


### 12 Channel, 21 mL Reservoir Plate

524092	21 mL	44	Rectangular	Pyramid	Natural	5/50
--------	-------	----	-------------	---------	---------	------


### 12 Channel, 290 mL Open format Reservoir Plate

524093	290 mL open format	44	Rectangular	Pyramid	Natural	5/50
--------	--------------------	----	-------------	---------	---------	------

Square Well 'U' Bottom    Rectangle Pyramid



Square Well Pyramid Bottom




## SAMPLE LIBRARY / 96 WELL CLUSTER TUBE 1.1 mL

Cat #	Description	Height (mm)	Well Shape	Well Bottom	Colour	Packing unit / Case
510083	8 Tubes Strip without Cap	46	Round	U	Natural	120
510084	Individual Tubes without Cap	46	Round	U	Natural	960
510085	Strip Cap for 8 Tubes strip	6	Round	-	Natural	120
510086	Racked 8 Tubes Strip, without Cap	49	Round	U	Natural	10 racks
510087	Racked 8 Tubes Strip without Cap Sterile	49	Round	U	Natural	10 racks



## SILICONE SEALING MAT FOR PLATES

Cat #	Description	Well Shape	Colour	Packing Unit / Case
510171	Sealing mat for 96 Well DeepWell Plate 0.8 ml / 1.2ml / Sample Library Tubes strip/ 390µl Micro Plates	Round	Natural	5 / 50
510172	Sealing mat for 96 Well DeepWell Plate 1.1ml / 2.0 ml / 0.5 ml	Round	Natural	5 / 50
510173	Sealing mat for 96 Well DeepWell Plate 2.2ml / 2.4 ml	Square	Natural	5 / 50
510174	Sealing mat for 48 Well DeepWell Plate 5 ml	Rectangle	Natural	5 / 50
510175	Sealing mat for 24 Well DeepWell Plate 10 ml	Square	Natural	5 / 50



## KING FISHER™ PLATES

### 96 Well KingFisher™ Style Plate


Cat #	Well Volume	Height (mm)	Well Shape	Well Bottom	Colour	Packing unit / Case
615020	2.0 ml	44	Square	Pyramid	Natural	5/50

### 96 Well KingFisher™ Style Micro Plate


615010	200 µl	14.5	Square	Pyramid	Natural	5/50
--------	--------	------	--------	---------	---------	------

### 96 Tip Comb





615030	-	43	Round	V	Natural	10/100
--------	---	----	-------	---	---------	--------



\*\* KingFisher is the registered Trade Mark of Thermo Fisher Group.



**TARSONS PRODUCTS LIMITED**  
902, Martin Burn Business Park,  
BP-3, Saltlake, Sector - V, Kolkata - 700091  
Tel: +91 033 - 35220300, e-mail: info@tarsons.com  
Website: www.tarsons.com

 TarsonsProductsPvtLtd  tarsons-products-ltd  tarsons.products  TarsonsProduct

