



cellpette®

Serological Pipettes



- Sterile & Non-cytotoxic
- Free of detectable Pyrogen, RNase & DNase
- Made of ultra-pure virgin polystyrene with aerosol resistant PET filter
- Calibrated to deliver (TD) with $\pm 2\%$ accuracy
- Negative graduations for extra capacity
- Manufactured in a Class 8 cleanroom
- Individually wrapped in thermoformed Paper / Plastic or Plastic / Plastic peel pack
- Ideal for cell culture, microbial culture and other research application



~Deliver with confidence

cellpette®



Cellpette, Sterile, Economy Pack

Individually wrapped in thermoformed paper/plastic peel pack

Code	Size(mL)	Colour code	Graduation(mL)	Tips	Tail	Packing
940010	1		1/100	Stretched	-	50 pc X10
940020	2		1/50	Stretched	-	50 pc X10
940030	5		1/10	Stretched	-	50 pc X 5
940040	10		1/10	Stretched	Welded	50 pc X 4
940140	10		1/10	Welded	Welded	50 pc X 4
940000	25		1/5	Welded	Welded	50 pc X 4
940001	50		1/2	Welded	Welded	25 pc X 4

Cellpette, Sterile, Economy Bulk Pack

Individually wrapped in thermoformed paper/plastic peel pack

Code	Size(mL)	Colour code	Graduation(mL)	Tips	Tail	Packing
94B001	1		1/100	Stretched	-	200 pc X 4
94B002	2		1/50	Stretched	-	150 pc X 4
94B005	5		1/10	Stretched	-	50 pc X 4
94B010	10		1/10	Stretched	Welded	50 pc X 4
94BW10	10		1/10	Welded	Welded	50 pc X 4
94B025	25		1/5	Welded	Welded	50 pc X 4
94B050	50		1/2	Welded	Welded	30 pc X 3

Cellpette, Sterile, Dispenser Box Pack

Individually wrapped in thermoformed paper/plastic peel pack

Code	Size(mL)	Colour code	Graduation(mL)	Tips	Tail	Packing
940011	1		1/100	Stretched	-	200 pc x 4
940021	2		1/50	Stretched	-	150 pc x 4
940031	5		1/10	Stretched	-	100 pc x 4
940041	10		1/10	Stretched	Welded	100 pc x 4
940141	10		1/10	Welded	Welded	100 pc x 4
940051	25		1/5	Welded	Welded	50 pc x 4
940061	50		1/2	Welded	Welded	40 pc x 4

Cellpette, Sterile, Dispenser Box Pack

Individually wrapped in thermoformed plastic/plastic peel pack

Code	Size(mL)	Colour code	Graduation(mL)	Tips	Tail	Packing
940012	1		1/100	Stretched	-	200 pc x 4
940022	2		1/50	Stretched	-	150 pc x 4
940032	5		1/10	Stretched	-	100 pc x 4
940042	10		1/10	Stretched	Welded	100 pc x 4
940142	10		1/10	Welded	Welded	100 pc x 4
940052	25		1/5	Welded	Welded	50 pc x 4
940062	50		1/2	Welded	Welded	40 pc x 4



For more information,
scan here



Serological Pipettes

Cellpette, Sterile, Bulk Pack

Code	Size(mL)	Colour code	Graduation(mL)	Tips	Tail	Packing
940013	1		1/100	Stretched	-	50 pc X 8
940023	2		1/50	Stretched	-	50 pc X 8
940033	5		1/10	Stretched	-	25 pc X 8
940043	10		1/10	Stretched	Welded	25 pc X 6
940143	10		1/10	Welded	Welded	25 pc X 6
940053	25		1/5	Welded	Welded	10 pc X 6
940063	50		1/2	Welded	Welded	10 pc X 5

Open End Cellpette, Sterile, Dispenser Box Pack

Individually wrapped in thermoformed paper/plastic peel pack

Ideal for use with high-viscosity liquids. Fast fill and release. Used with oils, cosmetics, paint, sludge, food, and cell culture applications.

Code	Size(mL)	Colour code	Graduation(mL)	Tips	Tail	Packing
940014	1		1/100	Open End	-	200 pc X 4
940024	2		1/50	Open End	-	150 pc X 4
940034	5		1/10	Open End	-	100 pc X 4
940044	10		1/10	Open End	Welded	100 pc X 4

Aspirating Cellpette, Sterile, Dispenser Box Pack

Individually wrapped in thermoformed paper/plastic peel pack. Without Aerosol - resistant PET filter

Economical, non-graduated, non-filtered pipettes. Ideal choice where precise measurement of transfer is not required. Safe alternative to glass aspirating pipettes.

Code	Size(mL)	Colour code	Graduation(mL)	Tips	Tail	Packing
940025	2	X	X	Stretched	--	150 pc X 4
940035	5	X	X	Stretched	-	100 pc X 4
940045	10	X	X	Stretched	Welded	100 pc X 4
940145	10	X	X	Welded	Welded	100 pc X 4

Bacteriological / Milk Cellpette, Sterile, Dispenser Box Pack

Available in single paper/plastic peel pack OR 25 pcs per plastic bag

Designed for use in the dairy market. Meet the Standard Methods for the Analysis of Dairy Products, APHA. Calibrated to deliver (TD) at +/-2%

Code	Size(mL)	Colour code	Graduation(mL)	Tips	Tail	Packing
940016	1.1 single pack		0.5,1.0,1.1	Stretched	-	200 pc X 4
940026	1.1 (25 / pack)		0.5,1.0,1.1	Stretched	-	25 pc X 32
940036	2.2 (Single pack)		0.5,1.0,2.0,2.1,2.2	Stretched	-	150 pc X 4
940046	2.2 (25 / pack)		0.5,1.0,2.0,2.1,2.2	Stretched	-	25 pc X 24




For more information,
scan here



cellpette®

Serological Pipettes



Cellpette, Sterile, Triple - Bagged

Individually wrapped in thermoformed paper/plastic peel pack are further packed in triple bag
Ideal for cell-based vaccine and drug production facilities, or another cleanroom laboratory for cell culture, microbial culture, and other research applications.

Code	Size(mL)	Colour code	Graduation(mL)	Tips	Tail	Packing
940017	1		1/100	Stretched	-	50 pc X 16
940027	2		1/50	Stretched	-	50 pc X 12
940037	5		1/10	Stretched	-	25 pc X 16
940047	10		1/10	Stretched	Welded	25 pc X 16
940147	10		1/10	Welded	Welded	25 pc X 16
940057	25		1/5	Welded	Welded	10 pc X 20
940067	50		1/2	Welded	Welded	10 pc X 15

Serological Pipette



*Cellpette is a registered Trademark of TARSONS PRODUCTS LIMITED

Cellpette_new_2025_02



TARSONS PRODUCTS LIMITED
902, Martin Burn Business Park,
BP-3, Saltlake, Sector - V, Kolkata - 700091
Tel: +91 033 - 35220300, e-mail: info@tarsons.com
Website: www.tarsons.com