





# Polypropylene Storage Plates




Microplates

Deep-Well Plates

KingFisher™ Plates

Sample Library Cluster Tubes

Sealing Mats



Designed for Safe Storage and High-throughput Screening




# Polypropylene Plates


Microplates | Deep-Well Plates | KingFisher™ Plates | Sample Library Cluster Tubes | Sealing Mats

- Plates are designed according to ANSI/SLAS format for high throughput screening, sample collection, compound storage, compound libraries, and combinatorial chemistry application.
- Compatible with most multichannel pipettes and automated Liquid-handling equipment.
- Made of virgin polypropylene, conforming USP class VI and compatible with solvents such as DMSO, phenol, and chloroform.
- Manufactured in ISO class 8 cleanroom facilities. Certified free of detectable DNase, RNase, and endotoxin.
- Alphanumeric well coding moulded in A-H and 1-12 axis.
- Raised rims on each well facilitate secure sealing. Compatible with adhesive seals or Silicone Sealing Mats.
- Temperature range for use: -80°C to 120°C.

Round-Well U-Bottom      Round-Well V-Bottom



Round-Well Flat-Bottom



## MICROPLATES

### 96-Well Standard Volume Microplates

Cat #	Well Volume	Height (mm)	Well Shape	Well Bottom	Colour	Packing unit/Case
510060	390 µL	14.5	Round	Flat	Natural	5/50
510060BK	390 µL	14.5	Round	Flat	Black	5/50
510063	390 µL	14.5	Round	Flat	White	5/50
510061	350 µL	14.5	Round	U	Natural	5/50
510061BK	350 µL	14.5	Round	U	Black	5/50
510064	350 µL	14.5	Round	U	White	5/50
510062	340 µL	14.5	Round	V	Natural	5/50
510062BK	340 µL	14.5	Round	V	Black	5/50
510065	340 µL	14.5	Round	V	White	5/50

### 96-Well Extended Volume Microplates

520071	500 µL	14.5	Round	V	Natural	5/50
520072	500 µL	14.5	Round	U	Natural	5/50

### 384-Well Microplates

520073	130 µL	14.5	Square	U	Natural	5/50
520074	130 µL	14.5	Square	Conical	Natural	5/50

## DEEP-WELL PLATES

### 96-Well Low-profile Deep-Well Plate for use with Magnetic Separator

Cat #	Well Volume	Height (mm)	Well Shape	Well Bottom	Colour	Packing unit/Case
510066	0.8 mL	31	Round	V	Natural	5/50

### 96-Well Low-profile Deep-Well Plates

510067	1.1 mL	27	Round	U	Natural	5/50
--------	--------	----	-------	---	---------	------

### 96-Well Standard-profile Deep-Well Plates

510069	1.2 mL	41.5	Round	U	Natural	5/50
510068	2.0 mL	45	Round	U	Natural	5/50
510077	2.2 mL	44	Square	Pyramid	Natural	5/50
510078	2.4 mL	44	Square	U	Natural	5/50

### 384-Well Low-profile Deep-Well Plates

510220	240 µL	22	Square	Conical	Natural	5/50
--------	--------	----	--------	---------	---------	------

~Technical specifications may be changed without prior notice.

# Polypropylene Plates

Microplates | Deep-Well Plates | KingFisher™ Plates | Sample Library Cluster Tubes | Sealing Mats

## HIGH VOLUME DEEP-WELL PLATES

### 48-Well, 5mL Deep-Well Plates

Cat #	Well Volume	Height (mm)	Well Shape	Well Bottom	Colour	Packing unit/Case
510080	5 mL	44	Rectangular	Pyramid	Natural	5/50
510080S	5 mL Sterile	44	Rectangular	Pyramid	Natural	5/50

### 24-Well, 10 mL Deep-Well Plate

510081	10 mL	44	Square	Pyramid	Natural	5/50
510081S	10 mL Sterile	44	Square	Pyramid	Natural	5/50
510082	10 mL	44	Square	U	Natural	5/50
510082S	10 mL Sterile	44	Square	U	Natural	5/50


### 12-Channel, 21mL Reservoir Plate

524092	21 mL	44	Rectangular	Pyramid	Natural	5/50
--------	-------	----	-------------	---------	---------	------


### 12-Channel, 290mL Open format Reservoir Plate

524093	290 mL open format	44	Rectangular	Pyramid	Natural	5/50
--------	--------------------	----	-------------	---------	---------	------

Square-Well U-Bottom      Square-Well Conical-Bottom



Rectangular Pyramid-Bottom



## SAMPLE LIBRARY / 96-WELL CLUSTER TUBE 1.1 mL

Cat #	Description	Height (mm)	Tube Shape	Tube Bottom	Colour	Packing unit/Case
510083	8-Tubes Strip without Cap	46	Round	U	Natural	120 Strips
510084	Individual Tubes without Cap	46	Round	U	Natural	960 Tubes
510085	8-Cap strip	6	Round	-	Natural	120 Strips
510086	Racked 8-Tubes Strip, without Cap	49	Round	U	Natural	10 Racks
510087	Racked 8-Tubes Strip without Cap Sterile	49	Round	U	Natural	10 Racks
510090	Empty rack with cover, without cap and Tube	52	-	-	Natural	10 Racks
510094	Racked individual Tubes, without Cap	49	Round	U	Natural	10 Racks



## SILICONE SEALING MAT FOR PLATES

Cat #	Description	Well Shape	Colour	Packing Unit/Case
510171	Sealing mat for 96-Well Deep-Well Plate 0.8 mL/1.2 mL/Sample Library Cluster Tubes/ 340/350/390 µL Microplates	Round	Natural	5/50
510172	Sealing mat for 96-Well Deep-Well Plate 1.1 mL/2.0 mL/ 500 µL	Round	Natural	5/50
510173	Sealing mat for 96-Well Deep-Well Plate 2.2 mL/2.4 mL	Square	Natural	5/50
510174	Sealing mat for 48-Well Deep-Well Plate 5 mL	Rectangle	Natural	5/50
510175	Sealing mat for 24-Well Deep-Well Plate 10 mL	Square	Natural	5/50



## KingFisher™ PLATES

### 96-Well KingFisher™ Style Plate

Cat #	Well Volume	Height (mm)	Well Shape	Well Bottom	Colour	Packing unit/Case
615020	2.0 mL	44	Square	Conical	Natural	5/50


### 96-Well KingFisher™ Style Microplate

615010	200 µL	14.5	Square	Conical	Natural	5/50
--------	--------	------	--------	---------	---------	------

### 96-Tip Comb

615030	-	43	Round	V	Natural	10/100
--------	---	----	-------	---	---------	--------

\*\* KingFisher is the registered Trademark of Thermo Fisher Group.



Tarsons\_Folder\_DeepWellPlates\_Feb\_2026



**TARSONS PRODUCTS LIMITED**

902, Martin Burn Business Park,  
BP-3, Saltlake, Sector - V, Kolkata - 700091  
Tel: +91 033 - 35220300, e-mail: info@tarsons.com  
Website: www.tarsons.com

 TarsonsProductsPvtLtd  tarsons-products-ltd  tarsons.products  TarsonsProduct

