

ACCUPENSE NEO

~Bottle Top Dispenser

~ Precision in Every Dispense

Key Features

- Positive displacement piston for accurate and reproducible dispensing
- Wide chemical compatibility for organic and inorganic solutions
- Leak-proof discharge valve with integrated safety bulb
- Fast and Precise volume adjustment
- 360° rotatable for ergonomic handling and clear bottle label visibility
- Plastic-Coated glass cylinder for enhanced user safety
- PFA sealing lip minimizes crystal build-up
- Easy recalibration for sustained accuracy
- Adjustable telescopic filling tube fits various bottle heights
- Fully autoclavable at 121 °C (2 bar) according to DIN EN 285

Material Used

- Cylinder- Borosilicate-Glass (Plastic-Coated)
- Piston Seal- PFA
- Wetted Parts- Borosilicate-Glass, PFA, PTFE, FEP, ETFE, PP, Al₂O₃, Ceramic, Platinum Iridium
- Thread Adapters- PP
- Filling Tube- Fluoroplastic

Specifications

Volume range (mL)	Inaccuracy	Imprecision
0.25 - 100.0	0.5%	≤ 0.2%

Ordering Information

Cat no.	Capacity (mL)	Graduation (mL)	Thread	Adapter
BTD-A002	0.25 - 2.5	0.05	GL45	GL 25-26, 28, 32 GL 32-33, 38, 40
BTD-A005	0.5 - 5.0	0.10	GL45	GL 25-26, 28, 32 GL 32-33, 38, 40
BTD-A010	1.0 - 10.0	0.20	GL45	GL 25-26, 28, 32 GL 32-33, 38, 40
BTD-A030	2.5 - 30.0	0.50	GL45	GL 25-26, 28, 32 GL 32-33, 38, 40
BTD-A060	5.0 - 60.0	1.00	GL45	GL 25-26, 28, 32 GL 32-33, 38, 40
BTD-A100	10.0 - 100.0	2.00	GL45	GL 25-26, 28, 32 GL 32-33, 38, 40

Scope of Delivery

Bottle Top Dispenser (GL 45 thread), 6 thread adapters, Adjustable telescopic filling tube, Quality certificate, Operating manual and Mounting tool

1
Year
Warranty



TARSONS PRODUCTS LIMITED

902, Martin Burn Business Park,
BP-3, Saltlake, Sector - V, Kolkata - 700091
Tel: +91 033 - 35220300, e-mail: info@tarsons.com
Website: www.tarsons.com



TarsonsProductsPvt.Ltd



tarsons-products-ltd



tarsons.products



TarsonsProduct

